

LOKI[®] PL

CANINE HYDROPHOBIC LENS

PRELOADED



- Monofocal
- For veterinary use



TECHNICAL SPECIFICATIONS

Lens type	For implantation in the canine capsular bag
Optic diameter	6.50 mm
Overall diameter	12.00 mm - 13.00 mm - 14.00 mm
Design	One-piece
Optic design	Monofocal Aspherical with negative spherical aberration to partly correct corneal spherical aberration Biconvex
Angulation	5°
Material	Hydrophobic CBK 1.8 from Cristalens
Dioptic power	+41.0D
Refractive index	1.54
Sterilization	Ethylene oxide
Injection system	Preloaded system
Recommended incision size	3.0 mm

