






PFI | **VET**



PFI S4 · PFI M4 · PFI X4

Aberration Balanced
Optimal Stability

Aberration Balanced



Type Single-piece Aspheric Intraocular Lens

Optic Aspheric - Coma Optimized

Material PFI VET Hydrophilic-Hydrophobic Copolymer
25% water content with UV Absorber

Optic diameter 6.5 mm

Overall diameter

- PFI S4 12.0 mm
- PFI M4 13.0 mm
- PFI X4 14.0 mm

Angulation 0°

Available powers +42.0D

PCO protection 360° Special Square Edge (patent pending)

Sterilization Steam

Implantation site Capsular bag