




CANABRAVA RING



Iris Expander Ring

<i>Brand name</i>	CANABRAVA RING
<i>Reference</i>	CANA - RING
<i>Description and Indications</i>	<p>The Iris Extender Ring is designed to provide stable and sustained pupillary dilation by mechanically extending the iris to facilitate the surgeon's field of vision and improve access to the posterior chamber during ophthalmic surgery.</p> <ul style="list-style-type: none">• 300° semicircular ring with a 60° opening made from a single piece.• Does not require an injector.• Easy to adapt and remove.• Easy to insert in pupils smaller than 4.5 mm and in shallow chambers.• Optimum visibility.• Single-use product <p>Indications:</p> <ul style="list-style-type: none">• Myosis or small pupils.• Iris coloboma.• Traumatic sectorial iridectomy.• Intraoperative floppy iris syndrome.• Prolapse of iris.• Tissue of the rigid iris due to a pseudoexfoliation.• Vitrectomy with small pupils.
<i>Material</i>	Medical grade Polymethyl Methacrylate (PMMA).
<i>Sterilization method</i>	Ethylene Oxide (ETO).
<i>Packaging</i>	Sterile, single unit.
<i>Shelf-life</i>	4.5 years

CE 0318



V-04-01-AAI-IMP REV.00