

ALOS

Palpebral Implants

Brand name

ALOS


Models

GOLD (TH) and PLATINUM (PTH)

Description and Indications

- Biocompatible materials
- Due to their curved design with rounded edges, the implants adapt to the curvature of the eyeball for a better esthetic result.
- Rounded edges to avoid the extrusion of the implant in the long-term.
- Suture holes (fixation channels gold references).
- Purity of 99,90%(Platinum) and 99,99%(Gold).
- Sterile.

Indication: Eyelid implants are designed for the correction of lagophthalmos.



Brand name **ALOS SIZER SET**

References **SET MEDIDOR PESAS TH**

Weight (g) 0.6-0.8-1.0-1.2-1.4-1.6

GOLD

References	TH 06	TH 08	TH 10	TH 12	TH 14	TH 16	TH 18
Weight (g)	0.6	0.8	1.0	1.2	1.4	1.6	1.8
Length (mm)	10.5	12.5	15.5	17.5	19.5	21.5	24.5
Width (mm)	5.7	5.7	5.7	5.7	5.7	5.7	5.7
Thickness (mm)	0.6	0.6	0.6	0.6	0.6	0.6	0.6
Sterilization method	Steam						
Shelf-life	3 years						

CE 2797

PLATINUM

References	PTH 06	PTH 08	PTH 10	PTH 12	PTH 14	PTH 16	PTH 18
Weight (g)	0.6	0.8	1.0	1.2	1.4	1.6	1.8
Length (mm)	10.2	12.4	15.0	17.0	18.8	21.0	24.0
Width (mm)	5.7	5.7	5.7	5.7	5.7	5.7	5.7
Thickness (mm)	0.6	0.6	0.6	0.6	0.6	0.6	0.6
Sterilization method	Steam						
Shelf-life	3 years						

CE 0318



V-04-01-ABQ-IMP REV.00