

AJL
vet

AJLVET VISCOELASTIC SOLUTIONS

DESCRIPTION AND INDICATIONS

Ophthalmic medical device designed for intraocular injection in surgeries in the anterior segment of the eye, such as cataract surgeries.

Maintains the anterior chamber's depth providing a high protection of the tissues and avoiding the risk of corneal endothelium trauma.

For veterinary use only.

HOW SUPPLIED

Sterile, in prefilled type I glass syringes of 2 ml with luer lock connection.

Each box contains five syringes and five cannulas.



MODELS & SPECIFICATIONS

HYALURONIC ACID (NaHA) SOLUTIONS

[ref. VV14] 1.4% 15,000-20,000 mPa·s (blue)
[ref. VV20] 2% 40,000-60,000 mPa·s (pink)
[ref. VV30] 3% 160,000-200,000 mPa·s (green)

Molecular weight: 2,400,000 Da

HPMC SOLUTIONS

[ref. VVCELL] 2% 2,000-40,000 mPa·s (clear)

Molecular weight: 92,000 Da